



CARE ABOUT MILLING

WWW.FREWITT.COM

# PMC-160

FREWITT, world leader in powder size reduction technologies, presents FlexMill, an innovative technology combining multiple milling processes in a single milling platform, dedicated to the pharmaceutical, chemical, battery and food product industries.

#### 2 PROCESSES - 2 MILLING HEADS - ONE SYSTEM

A unique, multi-process platform, the PMC-160 offers unequalled flexibility and extraordinary performance allows for significant savings.

Owing to a design enabling the use of 2 interchangeable milling heads on the same base module, FlexMill permits scale-up of 2 processes, 1 with or 1 without cryo – all within the same platform.

The PMC-160 is a high-performance pin mill designed for pulverizing hard or crystalline products to a particle size down to D90 of 10 µm. Ideal for processing small-scale laboratory batches of 50 g to 2 kg, or work in continuous mode up to 50 kg/h. The PMC-160 lends itself to be used for lab, pilot, or production scale needs. The device's compact design facilitates an easy integration and a quick set-up in any laboratory as well as an integration into production lines or isolators.

#### SUITABLE FOR THE FOLLOWING INDUSTRIES:











**Pharmaceutical** 

Battery

Chemical

Food

#### **MODULARITY**

A distinctly avant-garde, modular concept: FlexMill allows you to alternate, depending on your production requirements, between 2 milling heads for the micronization of low Tg products (glass transition temperature) and products resisting milling in ambient conditions.

#### **PERFORMANCE**

FlexMill proves its capabilities with the same degree of flexibility across the entire range of your processing needs, from small batch to production scale.

#### INNOVATION

Due to its extraordinary versatility and high performance, this mill can be used in many demanding milling processes, whether for achieving optimum particle size distribution, milling temperaturesensitive products, or potentially explosive substances.

#### **SAVINGS**

Enhanced return on investment and lower cost of ownership: FlexMill integrates two milling processes into a single platform. With its low level of consumption, the PMC-160 is the best solution for combining micronization and cost saving.

# **ADVANTAGES**

- Scalable to Frewitt Lab Scale / Production Scale Pin Mill systems
- Compact assembly for easy in-line integration
- Flexibility of modular arrangement
- Easy top-down cleaning (easy disassembly for easy cleaning)
- Reduced product contact surface areas
- Less energy consumption





**End fineness** Down to D 90 < 6 μm



Throughput Up to 50 kg/h



**PROCESS CHARACTERISTICS** 

Product temperature range -20°C to +60°C



**Rotor speed** Ultra high speed

# **EXECUTIONS**



# **Standard**

- Pharma execution GMP/FDA
- Fully automated system
- Integrated dosing system with anti-vault stirrer for entry funnel
- Integrated cooling system
- Pressure and temperature sensors
- Inline process
- Mobile cart with lifting device



# **Optional**

- ATEX
- WIP cleaning
- Controlled milling under inerted atmosphere
- Specific custom tooling
- 21 CFR part 11/ GAMP 5
- Flexible arrangement
- Automatic LN2 control

# **TECHNICAL CHARACTERISTICS**



# **Materials**

AISI 316/316L - Ra ≤ 0.8 µm AISI 304/304L - Ra ≤ 1.4 µm EPDM/PTFE/Silicone Interior Exterior In product contact, certified FDA





# **Dimensions & weight**

Length	1950 mm
Width	950 mm
Height	1800 mm
Weight	350 kg
Weight	350 kg



#### FREWITT HEADQUARTERS

# Frewitt manufacturer of machines Ltd.

Route du Coteau 7 1763 Granges-Paccot Switzerland

## Postal address

P.O. Box 1701 Fribourg Switzerland

P: +41 26 460 74 00 info@frewitt.com www.frewitt.com

# Frewitt Pharma Machinery (China) Co. Ltd.

1029 Shenfu Road, Xinzhuang Industrial Zone 201108 Shanghai China

P: +86 21 548 32950 info@frewitt-china.cn www.frewitt-china.cn

# Frewitt USA, Inc.

249 Homestead Road, Building 5, Unit 8 Hillsborough, NJ 08844 USA

P: +1 908 829 5245 info@frewitt-usa.com www.frewitt-usa.com